

DiSPIM – A Flexible Dual-View Light Sheet Microscope Platform



ASI's DiSPIM Team:

John Zemek	President
Gary Rondeau	Technical Director
Jon Daniels	DiSPIM Lead Engineer

NIH Collaborators and Inventors:

Hari Shroff	NIH/NIBIB Chief Scientist
Yicong Wu	NIBIB Staff Scientist
Abhishek Kumar	Postdoctoral Fellow


Micro-Manager Development:

Nico Stuurman	UCSF Vale Lab
Jon Daniels	ASI



DiSPIM = dual-view SPIM on inverted microscope



APPLIED SCIENTIFIC
INSTRUMENTATION



- Light sheet on inverted microscope
- Two (fixed) views → isotropic resolution
- Conventional sample mounting




Add-on for Conventional Inverted Microscope



- Existing mounts for:
 - Leica DMI-6000
 - Nikon TE-300, Ti
 - Olympus IX-71/81
 - Olympus IX-73/83
 - Zeiss Axio-Observer
 - Others can be designed

On ASI RAMM Frame

- Flexible Modular Inverted Microscope
- Environmental chamber option



Outline



APPLIED SCIENTIFIC
INSTRUMENTATION

- Introduction
- History
- Images
- Implementation details
- Micro-manager plugin
- Conclusion




APPLIED SCIENTIFIC
INSTRUMENTATION

History

- Hari Shroff & Yicong Wu (NIH)
 - Original iSPIM paper: Wu et.al., PNAS 108 (43) 17708-17713 (2011)
 - Original diSPIM paper: Wu et.al., Nature Biotechnology 31, 1032–1038 (2013)
- ASI contributions:
 - Modular microscope, including DiSPIM-specific parts
 - Fiber-coupled scanners
 - Controller-based synchronization
 - Micro-manager plugin


iSPIM Concept




- SPIM on inverted microscope
- Sample mounted on standard glass coverslip
- 30x faster than spinning disk for same SNR

Wu et. al, PNAS 108, 17708-17713 (2011)

DiSPIM Concept






Wu et. al, Nat. Biotech.
31, 1032-1038 (2013)



- iSPIM but use both objectives for imaging => dual orthogonal views => isotropic resolution after computationally combining

Isotropic Resolution




Joint Decon: A. York and Y. Wu

Wu et al. *Nat. Biotechnol.* 31, 1032-138 (2013), Kumar et al. *Nature Protocols* 9, 2555-2573 (2014), Ingaramo et al. *Chem Phys Chem* 15, 794-800 (2014)

Early diSPIM




APPLIED SCIENTIFIC
INSTRUMENTATION



Modern DiSPIM



APPLIED SCIENTIFIC
INSTRUMENTATION



1. SPIM mount
2. RAMM frame
3. Objective mount
4. Objective piezo
5. Bottom-side microscope
6. CDZ centering stage
7. SPIM LS-50 Z-drive
8. Camera mirror cubes
9. Excitation filter cubes
10. Camera tube lens
11. Scanner tube lens
12. Spacer
13. Light sheet scanners
14. XY stage (large MS2500)

Modern diSPIM



- Fits 90 x 75 cm air table
- Environmental control system
- Microscope alone costs <\$100k, additional \$75-150k for cameras, laser launch, epi light source, environmental control, computer, etc.



APPLIED SCIENTIFIC
INSTRUMENTATION

diSPIM Performance

- 3D volumes w/ isotropic resolution
 - Isotropic resolution requires registration/fusion
 - 330 nm resolution with Nikon 40X NA 0.8
- Acquisition rates up to 200 images per second (2-5 volumes per second)
- Achieve a $>10x$ reduction in photo-bleaching compared to con-focal methods.




APPLIED SCIENTIFIC
INSTRUMENTATION

Why DiSPIM?

- Sample mounting on coverslip/dish
- Isotropic resolution without rotating sample
 - Observe fast processes or moving objects
- Commercial but open and customizable
- Can be integrated with your favorite techniques or equipment


iSPIM vs. DiSPIM




Wu et. al, Nat. Biotechnol. 31, 1032-1038 (2013)

GFP histones in a live BV24 *C. elegans* embryo from the 4 cell stage up to hatching.
Sampled every minute at 50 planes per volume with 1 μm inter-plane spacing at 200 Hz.

DiSPIM vs. SDCM




- GFP-EB3 Microtubules in Live Human Umbilical Vein Endothelial Cells
- At same SNR and illumination, diSPIM enables collection of 3x more volumes, 3.2x more planes per volume and 7.6-fold less photobleaching



- C3D GFP-EB3 microtubule dynamics in human umbilical vein, endothelial cells
- Different thickness and cellular environments

Wu et. al, Nat. Biotechnol.
31, 1032-1038 (2013)





- Cultured human lung fibroblast cell, expressing GFP-tagged H-Ras, and demonstrating dynamic cellular movements

**Kumar et. al, *Nat. Protocols*
9, 2555–2573 (2014)**


Cancer Spheroids in Culture

View A

View B



Maximum projection vs. time, 2 views (not fused)
Courtesy Christian Conrad, Univ. of Heidelberg



- A549 cells infected for 16 hours with WSN PA-GFP and then imaged for 30 min with an entire cell volume captured every 2 seconds. Scale bar is 10 μm .


Lakdawala et al., PLoS Pathogen
10(3): e1003971 (2014).

Dual-color Imaging




APPLIED SCIENTIFIC
INSTRUMENTATION

- *C. elegans* embryo




Courtesy Abhishek Kumar
and Ryan Christiensen

Stage Scanning



Courtesy Abhishek Kumar



Scanned sheet



Stage scanning



APPLIED SCIENTIFIC
INSTRUMENTATION

ASI DiSPIM Implementation


- Built with modular hardware for evolving applications and configurations
- Modular control electronics for evolving needs
- Micro-Manager support w/ASI support
- Systems available from ASI directly or via multiple system integrators
- Many additional features/variations being developed by labs worldwide

ASI's Modular Microscope Hardware



APPLIED SCIENTIFIC
INSTRUMENTATION

- Adjustable cubes
- Tube lenses
- Camera mounts
- Illumination adapters
- Linear stages
- Filter wheels
- And many more LEGO-like components



Adjustable Beam-Splitter Cube



Four cameras





Transmitted light for bottom scope

Easy to build complex or custom microscope

Light Sheet Scanners

- Compact 2-axis scanner with fiber input
 - One axis makes sheet, other selects slice position
- C-mount device (19mm FN)
- 1kHz scanner bandwidth
- Anti-stripping option
- Tiger control card
- Many other applications



Tiger Controller

- Basic diSPIM has 10 controlled axes
 - XY stage
 - SPIM height
 - Lower Z
 - Two objective piezos
 - Two 2-axis scanners
 - Plus TTL triggers
- All controlled with ASI's Tiger modular controller
- Piezos and scanners controlled with either internal DACs or external voltage control



Hardware Synchronization

- Light sheet, piezos, cameras, and lasers must be tightly synchronized
- 2 approaches:
 - Generate synchronized control voltages
 - Labview code available from Shroff group
 - Third party software or write your own
 - Use synchronization within Tiger controller
 - Micro-Manger DiSPIM plug-in
 - Third party software or write your own

Micro-manager plugin features

- Acquisition Modes:
 - Single or double-sided
 - Synchronized slice and piezo
 - Fixed sheet
 - Stage scan
- Multi-Dimensional Acq.
 - Time points
 - Multi-position
 - Multi-color
- Supported cameras:
 - Andor Zyla
 - PCO Edge
 - Hamamatsu Flash 4
 - Others possible
- Supported lasers:
 - Lasers with dual port switch or passively split
 - Lasers on/off via TTL
 - Up to 4 colors

Navigation Tab

- Move all 10 axes
- Flexible manual control with joystick and wheels



ASI diSPIM Control

Navigation

- Setup Path A
- Setup Path B
- Acquisition
- Data Analysis
- Devices
- Settings
- Help

Joystick:

Left Wheel:

Right Wheel:

Path A: Beam Sheet

Path B: Beam Sheet

Change settings on tab activate


Camera:

Scan A
 Scan B
 Piezo A
 Piezo B

XY Stage, X axis: -2,819.6 μm	<input type="text" value="0"/>	-	<input type="text" value="10"/>	+	<input type="button" value="Go to 0"/>	<input type="button" value="Set 0"/>
XY Stage, Y axis: 690.52 μm	<input type="text" value="0"/>	-	<input type="text" value="10"/>	+	<input type="button" value="Go to 0"/>	<input type="button" value="Set 0"/>
Lower Z Drive: 5,792.94 ...	<input type="text" value="0"/>	-	<input type="text" value="100"/>	+	<input type="button" value="Go to 0"/>	<input type="button" value="Set 0"/>
Upper (SPIM) Z Drive: 20,000.0...	<input type="text" value="0"/>	-	<input type="text" value="100"/>	+	<input type="button" value="Go to 0"/>	<input type="button" value="Set 0"/>
Imaging Piezo A: 0 μm	<input type="text" value="0"/>	-	<input type="text" value="5"/>	+	<input type="button" value="Go to 0"/>	
Imaging Piezo B: 0 μm	<input type="text" value="0"/>	-	<input type="text" value="5"/>	+	<input type="button" value="Go to 0"/>	
Scanner A, sheet (X): 4 $^{\circ}$	<input type="text" value="0"/>	-	<input type="text" value="0.2"/>	+	<input type="button" value="Go to 0"/>	
Scanner A, slice (Y): 4 $^{\circ}$	<input type="text" value="0"/>	-	<input type="text" value="0.2"/>	+	<input type="button" value="Go to 0"/>	
Scanner B, sheet (X): 4 $^{\circ}$	<input type="text" value="0"/>	-	<input type="text" value="0.2"/>	+	<input type="button" value="Go to 0"/>	
Scanner B, slice (Y): 4 $^{\circ}$	<input type="text" value="0"/>	-	<input type="text" value="0.2"/>	+	<input type="button" value="Go to 0"/>	

Setup Tabs

- Useful for mechanical alignment
- Piezo/scanner cross-calibration



The screenshot displays the ASI diSPIM Control software interface. The window title is "ASI diSPIM Control". On the left, there is a navigation pane with the following tabs: Navigation, Setup Path A, Setup Path B, Acquisition, Data Analysis, Devices, Settings, and Help. The "Settings" tab is currently selected. Below the navigation pane, there are four status indicators: Scan A (black square), Scan B (red square), Piezo A (green square), and Piezo B (blue square).




The main control area is divided into several sections:

- Joystick:** A dropdown menu set to "XY Stage".
- Left Wheel:** A dropdown menu set to "Imaging Beam, slice (Y)".
- Right Wheel:** A dropdown menu set to "Imaging Piezo".
- Imaging side:** Checkboxes for "Beam" and "Sheet" (both checked) under "Imaging side", and "Beam" and "Sheet" (both unchecked) under "Epi side". A checkbox for "Change settings on tab activate" is checked.
- Camera:** A dropdown menu set to "LeftCam (Imaging)" and a "Live" button.
- Imaging center:** A field showing "0.0 μm " with "Go to" and "Set" buttons. A $\Delta = 5 \mu\text{m}$ field with up and down arrow buttons.
- Piezo = -0.978 μm + Slice * 108.947** with "Set calibration:" buttons for "2-point" and "Offset".
- Calibration Start Position:** "Slice position: 0 $^\circ$ " with a "0" input field and "Go to" button. "-0.5061 $^\circ$ " with a "Go to" button.
- Calibration End Position:** "0.5 $^\circ$ " with a "Go to" button. "Imaging piezo: 0 μm " with a "0" input field and "Set" button. "52.311 μm " with a "Set" button.
- Illum. piezo:** "0 μm " with a "0" input field, "Set home", and "Go home" buttons. A checked checkbox for "Go home on tab activate".
- Sheet width:** A slider control ranging from 0.0 to 8.0, with "-" and "+" buttons.
- Sheet offset:** A slider control ranging from -1.0 to 1.0, with "-" and "+" buttons, and a "Center" button.

Microscope Alignment




APPLIED SCIENTIFIC
INSTRUMENTATION




Co-aligning
objectives and
centering
cameras

Aligning
light sheets
hitting
cover-slip in
dye solution



Bottom Camera



MultiCam

Acquisition Tab

- Set parameters and acquire

The screenshot displays the ASI diSPIM Control software interface, specifically the Acquisition tab. The interface is organized into several sections:


- Navigation:** A vertical sidebar on the left contains buttons for Setup Path A, Setup Path B, Acquisition (selected), Data Analysis, Devices, Settings, and Help.
- Durations:** Shows acquisition parameters: Slice: 30.75 ms, Volume: 1,738 ms, Total: 1.738 s.
- Time points:** Includes a checkbox for Time points, a Number spinner set to 3, and an Interval [s] spinner set to 10.
- Multiple positions (XY):** Features a checkbox, an Edit position list... button, and a Post-move delay [ms] spinner set to 1.
- Channels:** Includes a checkbox, a Channel group dropdown set to Excitation Source, and a table for channel selection:

Use?	Preset	
<input checked="" type="checkbox"/>	488nm	+
<input checked="" type="checkbox"/>	561nm	-
- Volume Settings:** Contains various volume-related parameters: Number of sides (2), First side (A), Delay before side [ms] (500), Slices per volume (12), Slice step size [μm] (1), Minimize slice period (checked), Slice period [ms] (20), and Sample exposure [ms] (20.5). A Calculate slice timing button is also present.
- Data Saving Settings:** Includes a checkbox for Separate viewer / file for each time point, Hide viewer, Save while acquiring (checked), Directory root (\\ASI_Test\AcquisitionData), and Name prefix (test).
- SPIM mode:** A dropdown menu set to Synchronous piezo/slice scan (standard).
- Start! Button:** A green button with a right-pointing arrow, labeled "Start!".
- Status:** A message box indicating "No acquisition in progress."
- Use Navigation joystick settings:** A checked checkbox.

At the bottom left, there is a legend for the acquisition components: Scan A (black square), Scan B (red square), Piezo A (green square), and Piezo B (blue square).

Data Analysis Tab

- Export to MIPAV
- Open directly in Fiji Multi-View Reconstruction



Lots of options

- Many system configurations possible
- Multiple system integrators offering diSPIM
 - 3i is most advanced/proficient of integrators
 - Can also buy direct from ASI
- Multiple acquisition softwares
 - Including open source Micro-manager
 - 3i's Slidebook, others to come
- Multiple data analysis softwares
 - Including open source

Extensions


- Single-sided system with high-NA objective
- Photomanipulation from bottom
 - Optogenetics
 - Wounding
- 2-photon
- Imaging from bottom simultaneously
- Gated laser to avoid certain regions of sample

dispim.org wiki

← → ↻ dispim.org/start



Login



navigation

- ▼ docs
 - ▶ manual
 - comparison_of_spim_
 - comparison_with_confocal_
 - system_integrators
 - user_list
- ▼ hardware
 - add-ons
 - asi_manual
 - environmental_control
 - troubleshooting
- ▼ media
 - publications
 - video_tutorials
- ▼ software
 - fiji_mvr
 - labview
 - micro-manager
 - mipav_generatefusion
- registered users

search

toolbox

- What links here
- Recent changes
- Media Manager
- Site index
- Printable version
- Permanent link
- Cite this article

[article](#) [discussion](#) [show pagesource](#) [old revisions](#)

The Idea

The *diSPiM* is a flexible and easy-to-use implementation of Selective Plane Illumination Microscopy (*SPiM*) that allows for dual views (*d*) of the sample while mounted on an inverted (*i*) microscope. The *diSPiM* was co-developed by the NIH-based lab of Hari Shroff [and](#) Applied Scientific Instrumentation [\(ASI\)](#). *SPiM* is also referred to as light sheet fluorescence microscopy [or](#) LSFM because it uses a sheet or plane of light to illuminate the sample perpendicular to the imaging direction.

Table of Contents

- The Idea
- The Implementation
- The Benefits
- This Wiki

The Implementation

The *diSPiM* has two (usually symmetric) optical paths for light sheet imaging. Two objectives are placed at right angles above a sample mounted horizontally. A light sheet is created from one objective and then imaged using the other objective. A stack of images is collected by moving the light sheet through the sample. For a few users, the 3D information from a single view or stack is sufficient (*iSPiM*). For dual-view systems, the role of the two objectives is reversed to collect another stack from a perpendicular direction, and then the two datasets can be computationally merged to yield a 3D dataset with isotropic resolution (the usual problem of poor axial resolution is overcome by information from the other view).

The *diSPiM* "head" can be mounted on various inverted microscopes including ASI's RAMM frame (shown here). *diSPiM* systems can be obtained from various [system integrators](#). Various open-source and proprietary software packages are available for both data acquisition and data visualization. Most of the underlying microscope hardware is identical regardless of the system integrator and software used.

The choice of *diSPiM* objectives is limited because they must be co-focused without bumping into each other. The most commonly-used objectives for *diSPiM* are 40x water-dipping objectives with a NA of 0.8 (Nikon CFI Apo 40XW NIR).

Most often sCMOS cameras are used for *SPiM* imaging. There are working *diSPiM* systems with Hamamatsu Flash4, Andor Zyla, and PCO Edge cameras.

ASI makes a compact fiber-coupled 2D galvo or "scanner" which is an integral part of the system. It creates the light sheet by fast scanning in one axis and moves the sheet through the sample using the other axis. The output of the excitation laser (or laser launch) simply is fed into the scanner; it is helpful to have a 2x1 optical switch or dual-output laser launch so the excitation can all be steered to the scanner in the active light path.

For applications where [environmental control](#) is important, the *diSPiM* can easily be fitted with an incubator enclosure and appropriate equipment to keep samples alive and happy.

The bottom objective (the inverted microscope) most often has a lower-magnification objective and less expensive camera for locating the sample. Epi-illumination can easily be added.

The Benefits

Like other *SPiM* techniques, *diSPiM* illuminates only the focus plane and is thus ideal for imaging living cells and organisms because it minimizes photobleaching and phototoxicity effects.

Compared with confocal or spinning disk systems, there is ~2x improvement in axial resolution, >10x reduction in photobleaching, and speed comparable to spinning disk. See a more detailed [comparison with confocal](#).

A major advantage compared with most other *SPiM* implementations is that sample mounting is extremely simple, similar to an inverted microscope. Most commonly, specimens are placed on a 24x50mm coverslip which is held in a special chamber that holds the dipping media. See a more detailed [comparison with other SPiM techniques](#).





APPLIED SCIENTIFIC
INSTRUMENTATION

Thank you! Any questions?



APPLIED SCIENTIFIC
INSTRUMENTATION

Commercial SPIMs

<u>SPIM Type</u>	<u># views</u>	<u>Mounting</u>	<u>Software</u>	<u>Comments</u>
DiSPIM	2 fixed	Coverslip/dish with media	Open & various proprietary	Flexible configuration
Zeiss Z-1	Unlimited (rotation)	Capillary with agarose	1 proprietary	OpenSPIM is open alternative
Leica TCS SP8 DLS	1 fixed	Dish with media	1 proprietary	Combined with confocal
LaVision BioTec Ultramicroscope	1 fixed	??	1 proprietary	Mainly for large fixed samples